

# DIANA

## assay for inhibitor screening

### Overview and collaboration proposal

Institute of Organic Chemistry and Biochemistry AS CR, Prague, Czech Republic

---



**ÚOCHB** AV  
ČR  
**IOCB PRAGUE**



**IOCB TTO**  
[www.iocb-tto.cz](http://www.iocb-tto.cz)

## What is DIANA

- A **multi-well plate based assay** (with similar protocol to ELISA)
- Utilizing **detection probe consisting of small-molecule** active site ligand linked to reporter DNA oligo
- **Quantification by qPCR** with high sensitivity and broad dynamic range
- Can be implemented using standard equipment and be fully automatized - suitable for both academic labs and diagnostics / screening facilities

## Two major application markets

1. **DIANA for diagnostics**  
**Ultra-sensitive detection of active enzymes in range of clinical samples**  
(zeptomoles of target protein can be detected, >7-log dynamic range)
2. **DIANA for screening**  
**Screening for target enzyme inhibitors in drug discovery**  
(sensitive hit discovery, quantitative, cost-efficient)

# DIANA for screening: overview

## Application

- **Screening of compound libraries for new inhibitors / ligands** of target proteins
- Successfully tested and validated on multiple target proteins and protein families  
- assays for new targets straightforward to develop<sup>1</sup>

## Key advantages


- Extremely **high signal to noise ratio** ( $Z' > 0.9$ ;  $CV < 5\%$ )
- **Quantitative**: compound  $K_i$  directly measured from a single well, hits ranked by inhibition potency, easy to screen for specificity in protein family
- **Sensitive** : sensitive hit discovery, ultra-low false-positive and false-negative rate
- **Robust**: works with unpurified protein, no interference or non-specific binding
- **Cost-efficient**: very low compound consumption, allows for compound pooling leading to reduced screening costs and time

## Potential customers

- **Pharma and CROs screening facilities**: automated HTS drug discovery
- **Research institutions**: small-scale screens of in-house libraries against targets of interest (libraries of ~10k compounds can be screened in few 96-well plates)

1. Assuming at least one small-molecule inhibitor (even non-specific) exists and have suitable chemistry for modification.

# DIANA for screening: 5-step protocol easy to implement in most screening facilities




1. No temperature sensitive incubations. Flexibility in incubation times. Tolerates high DMSO contents.

# Key advantages: quantitative screening with broad dynamic range

## Quantification of affinity from a single read


- **Accurate determination of  $K_i$**  from measurement at single concentration of tested compound
- Screening **hits ranked by the inhibition** potency
- Possible to test **specificity against whole protein families** in a single run



Purified PSMA standard was titrated with PSMA inhibitor. Measured inhibition constants were stable over more than 5 orders of magnitude of inhibitor concentration

## Precise over broad concentration range

- Due to the qPCR detection, the same assay conditions can be used to precisely determine  $K_i$  of compounds **over ~7-log range** (e.g. 10pM - 100 $\mu$ M)
- DIANA measurements validated by orthogonal assays




$K_i$  of 41 different compounds with 7-log difference in PSMA affinity were determined using the same 100 $\mu$ M concentration. Measurements were validated by enzyme kinetics assays with serial dilution of inhibitors

# Key advantages: ultra-sensitive and robust screening assay

## Higher sensitivity of detection


- Extremely **high signal to noise ratio** ( $Z' > 0.9$ ,  $CV < 5\%$ )
- **Sensitive in hit discovery** ( $IC_{50} = K_i$ )
- **Ultra-low false-positive** and **false-negative rates**
- **Easy counter-screen setup** to exclude false negatives



Due to high sensitivity and low background, hits are identified in DIANA screening even in libraries, which failed to produce any hits in conventional screening techniques.

## Does not require recombinant protein

- Possible to screen inhibitors even with small amounts of targets endogenously present in body fluids – **no need for purified recombinant protein**
- Only **picograms of target protein** are sufficient to test the compound's activity<sup>1</sup>
- No fluorescent/colored compound interference




1  $\mu$ l of human sera containing approx. 5  $\mu$ g of endogenous PSMA was used for testing a compound's activity, with the same results as when using recombinant purified protein

1. Several orders less than what is usually required for ELISA or enzyme kinetics assays.

# Key advantages: compatibility with pooled compound screen offering high cost-efficiency

## Compound pooling setup...


- **Compounds pooled prior to initial screening** - only wells with strong hits subsequently analyzed per individual compound
- Enabled by **quantitative** nature of the assay, ultra-high sensitivity (screening at low concentrations) and ultralow false positive and false negative rate
- **Pilot study:** strongest hits in a library of >2000 compounds identified using only two 96-well plates<sup>1</sup>



- **In progress:** screening of ~25,000 compounds in a single 384-well plate with 80 compounds per well

## ... leading to highly cost-efficient screening

- Significant **reduction of time required** to complete the screening project
- Unparalleled **low screening costs** (also considering low amount of protein and probe required)



- **Enhanced screening capabilities** even for small-scale academic and R&D labs:
  - Typical academic research library of **2000 compounds** can be screened for inhibition potency to target using a kit containing only **two 96-well plates**<sup>1</sup>
  - Small screening facility with a single 384-well plate qPCR machine is able to screen up to **~0.5 million compounds per day**<sup>2</sup>

1. First 96 wells-plate: initial screening run with 24 compounds per well. Second 96 wells-plate: subsequent analyses of 4 strongest pools for  $K_i$  of individual compounds  
2. 16 plates per day (30 min qPCR) x 320 wells per plate x 80 compounds per well = 16 x 25600 = 409,600 compounds

# Catalogue of targets: expanding range of assays for commercially relevant targets

Status	Target	Relevant therapeutic area	Screening projects
Assay ready for screening	Prostate specific membrane antigen (PSMA)	<ul style="list-style-type: none"> <li>Oncology (Prostate cancer)</li> <li>Neurology (CNS)</li> </ul>	<ul style="list-style-type: none"> <li>Search for new scaffolds with desired pharmacokinetics</li> </ul>
	Carbonic anhydrase family (CA-II, VII, IX, XII)	<ul style="list-style-type: none"> <li>Oncology</li> </ul>	<ul style="list-style-type: none"> <li>Search for new specific inhibitors with desired pharmacokinetics</li> </ul>
	Fibroblast activating protein (FAP)	<ul style="list-style-type: none"> <li>Oncology</li> <li>Metabolic disorders</li> </ul>	<ul style="list-style-type: none"> <li>Screen for new scaffolds</li> </ul>
	FcγRI receptor, Glutamate carboxypeptidase III, Influenza neuraminidase, MTH1, ...		
Priority pipeline	Methyl transferases (e.g. EZH2)	<ul style="list-style-type: none"> <li>Oncology</li> <li>Epigenetics</li> </ul>	<ul style="list-style-type: none"> <li>Screen for specific inhibitors (no off-targets)</li> </ul>
	Hydroxysteroid dehydrogenase	<ul style="list-style-type: none"> <li>Woman health (Endometriosis)</li> <li>Oncology (Breast cancer)</li> </ul>	<ul style="list-style-type: none"> <li>Search for new inhibitor scaffolds</li> </ul>
	Kinases	<ul style="list-style-type: none"> <li>Multiple therapeutic areas</li> </ul>	<ul style="list-style-type: none"> <li>Pan-kinase probe allowing screening of inhibitors for any kinase target</li> </ul>
	Insulin receptor family	<ul style="list-style-type: none"> <li>Hormonal disorders</li> <li>Oncology</li> </ul>	<ul style="list-style-type: none"> <li>Screen for specific ligands</li> </ul>
On demand targets	Assay can be promptly developed on demand to majority of relevant protein targets <sup>1</sup>		

1. Assuming at least one small-molecule inhibitor (even non-specific) exists and have suitable chemistry for modification.



# Our proposition: looking for partnerships to introduce DIANA to drug discovery market

	Our proposition	Potential partners
1 Assay development and licensing for screening projects	<ul style="list-style-type: none"><li>• We <b>support implementation</b> of assay at <b>partner's site</b> and its adjustment to the screening project needs</li><li>• We <b>synthesize the detection probe</b> and provide other reagents required</li><li>• We can <b>develop DIANA based screening assay for any new target</b> of partner's interest</li></ul>	<ul style="list-style-type: none"><li>• Pharma screening center<sup>1</sup></li><li>• CRO<sup>2</sup></li><li>• Academic screening facility</li></ul>
2 HTS screening services	<ul style="list-style-type: none"><li>• We <b>screen partner's compound libraries</b> for hits to DIANA compatible targets and return quantification of their inhibition potency</li><li>• We can run <b>mid- to high-throughput screening projects</b> at our facility (up to 10<sup>6</sup> compounds in a pooled compound setup)</li></ul>	<ul style="list-style-type: none"><li>• Pharma</li><li>• Academic research groups<sup>3</sup></li><li>• Screening facility network</li></ul>
3 Distribution of 'screening kits'	<ul style="list-style-type: none"><li>• We <b>co-develop ready-made screening kits</b> and optimized instrumental setup – <i>see details in Appendix</i></li><li>• We <b>distribute via partner's product catalogue</b> to research laboratories and screening facilities for application in small-scale screening experiments<sup>3</sup></li></ul>	<ul style="list-style-type: none"><li>• Lab technology / chemical distributor</li><li>• Instrument manufacturer</li></ul>

# Contacts for enquiries & other resources

## Contacts

### Václav Navrátil – Head of R&D

- [vaclav.navratil@uochb.cas.cz](mailto:vaclav.navratil@uochb.cas.cz)
- Methodology questions, assay development requirements

### Jaromír Zahrádka & Martin Dienstbier – TTO

- [zahradka@iocb-tto.cz](mailto:zahradka@iocb-tto.cz), [dienstbier@iocb-tto.cz](mailto:dienstbier@iocb-tto.cz)
- Business questions, collaboration deals

### Other contacts

- Martin Fusek, Jan Konvalinka

## Further resources

### Research paper

- Navratil et al. (2017): DIANA for sensitive and selective enzyme detection and inhibitor screening. *Nucleic Acids research*
- Screening methodology & Clinical diagnostics papers. *In preparation*

### Introductory video

- [www.youtube.com/watch?v=hrh82euICfU](http://www.youtube.com/watch?v=hrh82euICfU)

Institute of Organic Chemistry and Biochemistry, Czech Academy of Sciences.  
Flemingovo namesti 2, 16610, Prague 6, Czech republic



ÚOCHB <sup>AV</sup><sub>ČR</sub>  
IOCB PRAGUE



IOCB TTO  
www.iocb-tto.cz

# Appendix: DIANA-based kits for small-scale screening can be developed for range of targets

## Application

- **Screening of small- to mid-size compound libraries** for new inhibitors / ligands of target protein (research laboratories, academic and pharma screening facilities)

## Kit content

- Target protein **coated 96- or 384-well plate** (or anti-GST or NTA coated + tagged target protein)
- Target specific **detection probe**
- Reconstitution buffers (for protein and probe)
- 10x assay buffers (protein and probe dilution, wash)

## Protocol and equipment required

- **Simple 3-step protocol:** incubation with compound and with probe in one step, washing, qPCR detection
- Requires **standard lab equipment** and **multiwell-plate qPCR system** (automated liquid handler and plate washer are optional)

