


BTES datasheet

Product #: 1

Product Name: BTES

Alternative names: 3-(4-((bis((1-(tert-butyl)-1H-1,2,3-triazol-4-yl)methyl)amino)methyl)-1H-1,2,3-triazol-1-yl)propane-1-sulfonic acid

Chemical structure:



Purity: 99%+ (determined by UPLC-MS)

Molecular weight: 494.62

Chemical composition: C₂₀H₃₄N₁₀O₃S

CAS code: N/A

Appearance: White to off-white powder

Solubility: Water, DMSO, DMF, NMP, MeOH, Acetone, CAN

Storage: 4 °C