


BTTP datasheet

Product #: 2

Product Name: BTTP

Alternative names: 3-(4-((bis((1-(tert-butyl)-1H-1,2,3-triazol-4-yl)methyl)amino)methyl)-1H-1,2,3-triazol-1-yl)propan-1-ol

Chemical structure:



Purity: 99%+ (determined by UPLC-MS)

Molecular weight: 430.56

Chemical composition: C₂₀H₃₄N₁₀O

CAS code: N/A

Appearance: White powder

Solubility: Water, DMSO, DMF, NMP, MeOH, Acetone, CAN, EtOAc

Storage: RT